



# Tutivia™ Test Workflow


1




2



3




4



5



6



1. The Tutivia test is ordered and specimen is collected to send to Verici Dx.  
2. Verici Dx laboratory isolates RNA, prepares library and conducts sequencing.

3. The sequencing data is uploaded and QC is assessed before files enter pipelines.  
4. Verici Dx proprietary pipelines and algorithm process data to produce test result.

5. Laboratory Director reviews assay and patient QC to approve results release.  
6. Verici Dx laboratory provides Tutivia results report to ordering clinician.